

The logo for LINC (Lung and Inflammation Network for Critical Care) features the letters 'LINC' in a white, sans-serif font. The letters are positioned over a stylized graphic of three curved, overlapping brushstrokes in dark blue, red, and yellow. The background of the slide is a light blue gradient with a large, faint, light blue brushstroke that curves across the upper left and middle sections.

LINC

Cardiac- and Respiratory-Induced Motion of Renal Arteries and Stents in Snorkel Endovascular Aneurysm Sealing

Christopher Cheng PhD

John Kim, Ga-Young Suh PhD, Andrew Holden MD

Disclosure

Speaker name: Christopher Cheng

.....

I have the following potential conflicts of interest to report:


- Consulting
- Employment in industry
- Stockholder of a healthcare company
- Owner of a healthcare company
- Other(s)

- I do not have any potential conflict of interest




Objective and Imaging

Quantify cardiac- and respiratory-induced renal artery and stent deformation in AAA patients treated with Sn-EVAS

10 patients underwent cardiac-resolved CTA during inspiration and expiration breath-holds





Modeling and Quantification



Results – Renal Snorkel Stent

(Values in Mean±SD)	Cardiac Deformation (Diastole to Systole)		Respiratory Deformation (Inspiration to Expiration)	
	Right	Left	Right	Left
Δ Arclength (%)	0.1±0.8	-0.4±1.2	0.3±1.8	-0.2±1.4
Δ Avg C (cm ⁻¹)	0.02±0.03	0.01±0.03	0.01±0.04	0.00±0.02
Δ Peak C (cm ⁻¹)	0.04±0.07	0.02±0.07	0.03±0.14	0.02±0.09
Max Δ C (cm ⁻¹)	0.13±0.09	0.15±0.11	0.20±0.10	0.18±0.09



- Significant cardiac-induced bending of the snorkel stent
- Significant respiratory-induced bending of the snorkel stent
- No significant difference between cardiac vs. respiratory deformation



Results – Renal Artery



(Values in Mean±SD)	Cardiac Deformation (Diastole to Systole)		Respiratory Deformation (Inspiration to Expiration)	
	Right	Left	Right	Left
ΔBranch Angle (°)	0±1	0±1	1±4	-3±4
ΔEnd-Stent Angle (°)	-1±2	0±2	3±5	10±7
ΔAvg C (cm ⁻¹)	-0.02±0.02	0.03±0.02	-0.01±0.11	0.06±0.10
Max ΔC (cm ⁻¹)	0.21±0.18	0.28±0.17	0.42±0.32	0.57±0.25

- End-stent angle change greater for respiration vs. cardiac
- Maximum curvature change greater for respiration vs. cardiac
- End-stent angle change greater for LRA vs. RRA for respiratory
- Larger variation in respiratory-induced deformations vs. cardiac



Conclusions

- Balloon-expanded snorkel renal stents used for Sn-EVAS undergo similar magnitude of bending deformations due to cardiac pulsatility and respiration
- The renal artery itself bends significantly more at the distal end of the stent with respiration as compared to cardiac by a factor of 2x.
- The left renal artery experiences greater bending while the right renal is relatively stabilized by the IVC.
- Renal arteries exhibit larger population variance for respiratory vs. cardiac deformations
- Cardiac and respiratory influences may challenge the mechanical durability of snorkel stents of Sn-EVAS similarly, however, respiration may be the primary culprit for tissue irritation, potentially increasing risk for end-stent thrombosis, especially in the left renal artery.





Thank You

The logo for LINC (Lung and Inflammation Network for Critical Care) features the word "LINC" in white capital letters. The letters are positioned over a stylized graphic of three curved, overlapping brushstrokes in dark blue, red, and yellow. The background of the slide is light blue with a large, faint, light blue brushstroke graphic that curves across the upper left and middle sections.

LINC

Cardiac- and Respiratory-Induced Motion of Renal Arteries and Stents in Snorkel Endovascular Aneurysm Sealing

Christopher Cheng PhD

John Kim, Ga-Young Suh PhD, Andrew Holden MD